



-25°C Medical Freezer

Instruction Manual

NOTE: Reasonable product improvements and technical updates without notice.

Dear Customer:

Congratulations on your purchase of this product. Like all of our products, this product has also been developed on the basis of the latest technology and produced using reliable and modern electrical components.

Please take a few minutes before starting operation of the appliance and read the following operating instructions. Save this manual for further reference.

Many thanks.





CAUTION!

All internal adjustments and maintenance must be performed by qualified service personnel.

This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.

Children shall not play with the appliance.

Cleaning and user maintenance shall not be made by children without supervision.



Environmental Protection

Discarded electric appliances are recyclable and should not be discarded in the domestic waste! Please actively support us in conserving resources and protecting the environment by returning this appliance to the collection centers (if available).



Important operating and/or maintenance instructions. Read the accompanying text carefully.



Potential electrical hazards. Only qualified persons should perform procedures associated with this symbol.



Hot surface(s) present which may cause burns to unprotected skin or to materials which may be damaged by elevated temperatures.



Extreme temperature hazards, hot or cold. Use special handling equipment or wear special, protective clothing.

Always use the proper protective equipment (clothing, gloves, etc.)

Always use dissipate extreme cold or heat and wear protective clothing.

Always follow good hygiene practices.

Each individual is responsible for his or her own safety.

Table of Contents

1.	Installation	2
1.1	Safety Information	2
1.2	Assembly	2
1.3	Positioning the Medical Freezer	3
1.4	Ventilation	4
1.5	Power Supply.....	4
2.	Product Description	5
3.	Operation And Display	5
3.1	Access to Power	5
3.2	Turn On/Off the Light.....	5
3.3	Indicator Light Status Instruction	6
3.4	Keys Instruction.....	6
3.5	Fault Code.....	7
3.6	Keys Operation - Parameter Adjustment.....	7
4.	Cleaning	8
4.1	Out of Service	8
5.	Maintenance.....	9
6.	Disposal.....	9

If the user's manual is lost or damaged, you can through our website to download or contact us by E-mail.


1. Installation

1.1 Safety Information

All electrical work must be executed by suitably qualified persons. When using any electrical appliance, safety precautions must always be observed.

- ◆ Do not use this appliance for other than its intended use. This product is intended for use as a Medical Freezer. Consult your supplier for alternatives.
- ◆ Do not cover the front grille or block the rear air entry by placing object up against the cabinet.
- ◆ Please close attention, when used by or near children, infirm persons.
- ◆ Ensure adequate ventilation.
- ◆ Do not probe any openings.
- ◆ Do not touch any moving parts or hot surfaces.
- ◆ Regulations require that all electrical work be carried out by authorized persons. For your own safety and that of others please ensure this is done.
- ◆ Do not overload the power supply. Always ensure on the rating label. The rating label is located. If have any problem, please consult a qualified electrician. Caution: Never use extension leads or multiple adaptors.
- ◆ If the mains supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified and skilled persons in order to avoid hazard.
- ◆ Disconnect the mains power supply before attempting any cleaning, removal of any covers, or maintenance work.
- ◆ Do not store explosive substances, such as aerosol cans with flammable propellant.
- ◆ Warning: Do not use mechanical devices or other means to accelerate the defrosting process. Defrosting is performed manually.
- ◆ Warning: Do not damage the refrigeration circuit.
- ◆ Warning: Do not use electrical appliances inside the storage compartment of this appliance.

1.2 Assembly

- ◆ Remove the appliance from the packaging and peel off any protective film from all surfaces.
 - ◆ Fit the shelf support clips onto the shelf support strips, making sure that the clips for each shelf are of the same height and that they are securely engaged on the support strip.
 - ◆ Unwrapping each shelf carefully, to prevent damaging their protective coating or surface finish.
- 

- ◆ Position the shelves onto the clips.



Please dispose of all waste packaging appropriately and in accordance with all local recycling laws.



WARNING

Do not overload the power supply. Always ensure that the power supply is adequate for the electrical requirements specified on the rating label and never use extension leads or multiple adaptors.

1.3 Positioning the Medical Freezer

- ◆ Please position the cabinet carefully being careful not to bump it.
- ◆ The power supply flex exits the Medical Freezer at the bottom right hand side (as viewed from the front). For ease of access; the flex should be retrieved prior to positioning the Medical Freezer in its final location.
- ◆ The maximum recommended ambient temperature for this Medical Freezer is 38°C
- ◆ Avoid positioning the appliance in direct sunlight or damp areas.
- ◆ Allow adequate space for the door(s) to fully open.
- ◆ Ensure that this Medical Freezer is positioned on a level surface, so as to allow the door(s) to shut and seal correctly. If necessary, adjust the screw feet of the appliance to make it level.

Note: If the unit has not been stored or moved in an upright position, let it stand upright for approximately 12 hours before operation.

1.4 Ventilation

- ◆ All models must have clear and unobstructed ventilation from the entire surface area of the front grille.
- ◆ The 2 and 3 door models may be built into an enclosure without the need for ventilation other than at the front.
- ◆ The single door of Medical Freezer must have a minimum ventilation space of 20mm at sides and rear.



- ◆ WARNING: Failure to provide the minimum ventilation space will harm the performance of your Medical Freezer could cause fire and will invalidate the warranty.




- ◆ WARNING: Keep clear of obstruction all ventilation openings in the appliance enclosure or in the structure for building in.

1.5 Power Supply

The Medical Freezer uses a power plug which must be well grounded. Never ground the Medical Freezer through a gas pipe, tap water pipe OR the grounded wire of single-phase power supply etc. Never share same circuit with other appliances. Before plugging in your Medical Freezer for the first time make sure the power cable is free from damage and check that your voltage corresponds to that on the serial plate.

Once plugged in the display will show the temperature inside the Medical Freezer and the compressor will run late on (the indicator light of cooling will be on at the same time). This temperature is the actual temperature recorded by the probes inside the Medical Freezer and as the fridge cools to the set point (factory set at 4°C/39°F) so the number on the display will descend accordingly. If power is cut off, wait for at least 5 minutes before plugging the unit in again to avoid damaging the compressor.

2. Product Description



3. Operation And Display

This product is designed to operate at a recommended ambient temperature not more than 38°C.

3.1 Access to Power





- ◆ Close the door of the appliance.

3.2 Turn On/Off the Light







- ◆ The lamp is controlled by the lamp switch.
- ◆ If you want to turn on the light while you do not open the door, you can turn on the switch.



3.3 Indicator Light Status Instruction

Indicator Light	Symbol	Status	Meaning
Setting	set	ON	Parameter setting
		OFF	Status of temperature measuring and controlling
Refrigeration		ON	Refrigeration work
		OFF	Refrigeration stop
		FLASH	Refrigeration time delay
Defrost		ON	Defrost work
		OFF	Defrost stop
Fan		ON	Fan work
		OFF	Fan stop
Defrost dripping	drip	ON	Start defrost dripping
		OFF	Stop defrost dripping
Door switch		ON	Cabinet door open
		OFF	Cabinet door close











3.4 Keys Instruction


Keys	Function	Button Action
Set	Enter the status of parameter setting;	Pressing the keys for 3s
	Switch between menu and parameter;	Press the response
	Adjust menu and parameters;	Press the response
	Open/close light (only valid for the model with light control)	Press the response
	Upload the data to copy card	pressing the keys for 3s
	Adjust menu and parameters	Press the response
	Download the copy card	pressing the keys for 3s
	View evaporator sensor temperature	Press the response
	Exit from parameter setting;	Press the response
	Press 3s to forced switch between refrigeration, defrost/defrost delay, defrost dripping	pressing the keys for 3s
Set + 	Check the maximum alarm temperature	Press the response
Set + 	Check the minimum alarm temperature	Press the response
Set + 	Clear alarm temperature record	Pressing the keys for 3s

3.5 Fault Code

Alarm code	Alarm reason
E1	control sensor failure
E3	Condenser sensor failure
E4	Cabinet temperature sensor failure
E00	Clear alarm temperature record
blinking Display	Power outage
cH	Condenser high temperature alarm
rH	Cabinet high temperature alarm
rL	Cabinet low temperature alarm
Er	Copy card programming failure
EP	Inconsistent data between copy card and controller, programming failure

3.6 Keys Operation - Parameter Adjustment

- 1) In the status of temperature measuring and controlling, press **Set** key for three seconds to enter user menu, it displays the code **St**, then press **Set** key again, display the value of St. It could be modified by pressing the key  or .
- 2) When it displays the code **St**, press the key , display the code **Po**, then press **Set** key, display **00**, at this time, press  or  to input the password of administrator menu.
- 3) Press **Set** key again to confirm the password input, and the controller will automatically verify the correctness of password. When it passes, it could select parameter items St、Po、C1、C2.....U1 (that is, any parameter items both in the administrator menu and user manuals) by pressing the key  or . Or else, only the parameters items **St** and **Po** available, others could not be displayed.
- 4) When the parameter item is selected, press **Set** key to enter to the setting of the current item, press  or  to modify the value, and then press Set key to return to the menu.
- 5) Under the status of parameter setting, press  key or no key operation within 30s, it will exit from parameter setting and automatically save the current parameter value.

Note: The password input of administrator menu only is valid for single entering. After exit from the parameter setting by pressing , it needs to input the correct password again for next parameter adjustment.

4. Cleaning



- ◆ **WARNING:** You must disconnect the plug before cleaning or maintenance.
- ◆ Do not use abrasive cleaning agents which may scratch and harm the delicate interior finishes of your appliance.
- ◆ Always wipe dry after cleaning.
- ◆ An agent or qualified technician must carry out repair if required.
- ◆ When cleaning the storage box, please use a soft cloth or sponge dipped in warm water (neutral detergent can be added) to clean the storage box.
- ◆ It is strictly forbidden to spray and wash with water to avoid affecting the electrical insulation performance.
- ◆ It is strictly prohibited to use hard brush, wire brush, detergent, gasoline, banana water, acetone and other organic solvents, hot water or acid, alkali, etc. to clean the storage box.

4.1 Out of Service

- ◆ If the storage box is out of service for a long time, the power supply should be cut off and cleaned according to the above method, and the door should be opened to dry.
- ◆ In general, it is not necessary to stop using the storage box for a long time in order to avoid affecting the service life.

5. Maintenance

If your appliance develops a fault, please check the following table before making a call to the helpline.

Fault	Possible Cause	Action
The appliance is not working	The unit is plugged in correctly	Check the unit is plugged in correctly
	Plug or lead is damaged	Call Our agent or qualified technician
	Power supply	Check power supply
	Internal wiring fault	Call Our agent or qualified technician
The appliance turns on, but the temperature is too high or too low	Filter or condenser blocked with dust	Clean filter or condenser
	Doors are not shut properly	Check door is shut and seals are not damaged
	Appliance is located near a heat source or air flow to the condenser is being interrupted	Move the appliance to a more suitable location
	Ambient temperature is too high	Increase ventilation or move appliance to a Medical Freezer position
	Insufficient airflow to the fans	Remove any blockages to the fans
	Appliance is overloaded	Reduce the amount stored in the appliance
	Factory default parameters adjusted	Call Our agent or qualified technician
The LED lights not working	Led light short Leaded damage	Call Our agent or qualified technician
The appliance appears BT alarm	See description below	See description below
The appliance is unusually loud	The appliance is touching a neighboring object	Check installation position and change if necessary
	The appliance has not been installed in a level or stable position	Check installation position and change if necessary

6. Disposal

If in any doubt, please consult your Local Authority, contact us or contact the reseller from whom you purchased the appliance

Certificate of Quality

Checker: