



-60°C Low Temperature Freezer Instruction Manual

NOTE: Reasonable product improvements and technical updates without notice.

Dear Customer:

Thanks for selecting DW-60W Series Low Temperature Freezer. We would provide the best product and service to you. Failure to read, understand and follow the instructions in this manual may result in damage to the unit, injury to operating personnel, and poor equipment performance.



CAUTION

- ◆ All internal adjustments and maintenance must be performed by qualified service personnel.
- ◆ This appliance can be used by persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- ◆ Children shall not play with the appliance.
- ◆ Cleaning and user maintenance shall not be made by children without supervision.



WARNING

1. As the freezer contains a refrigerant which are flammable, it is essential to ensure that the freezer is far away from any heat source.
2. Keep clear of all obstacles surrounding the freezer, and place the freezer in a well-ventilated area.
3. Do not store any material that is inflammable, explosive and corrosive inside of the freezer
4. Do not use mechanical devices or other means to accelerate the defrosting process, other than those recommended by the manufacturer.
5. Do not damage the refrigeration circuit by yourselves.
6. Do not put electrical devices inside the compartments of the freezer, unless they are of the type recommended by the manufacturer.
7. This freezer is not intended for use by young children or infirm person unless they have been adequately supervised by a responsible person to ensure they can use the freezer safely. Young children should be supervised to ensure that they do not play with the freezer.

8. Keep always the keys in a separate place and out of reach of children.
9. Never place containers with liquids, especially fizzy drinks, in the freezer. The liquid may cause an explosion during the freezing process.
10. Never touch the refrigerating surfaces with bare hands. Always wear gloves to open the freezer and load/unload the frozen items to prevent cold injury.



BEFORE USE


1. *Before operating the freezer, please read the following instructions carefully as they contain important information on safety, installation, operation and maintenance. Please keep the instructions at your hands for future possible usage.*
2. *While unpacking the freezer, please check to ensure the freezer has not been damaged and have all accessories packed well.*
3. *Transportation damage should be reported to the local distributor before the freezer is put into use.*
4. *Remove all the packaging material, and position the freezer in a stable, level and dry place.*
5. *Before the 1st time use, the freezer should be stood for at least 12 hours.*

Table of Contents

1. Installation and Start-up.....	6
1.1 Placement.....	6
1.2 Room Temperature.....	6
1.3 Ventilation.....	6
2. Electrical Connection	6
3. Temperature Controller	7
3.1 Description	7
3.2 List of Parameters	7
3.3 Parameter programming	7
3.4 Access to all code protected parameters	8
3.5 Default Working.....	8
3.6 Led Indication and Display Messages.....	8
4. Defrosting	8
5. Cleaning and Maintenance.....	9
5.1 Cleaning	9
5.2 Energy Saving	10
6. Trouble Shooting	10
7. Safeguarding on the Environment	11
7.1 Packing.....	11
7.2 Product.....	11
7.3 Disposal.....	11

The operating instructions should always be available and accessible for the person who use the freezer!

The product complies with CE marking regulation, Directives and relative standards.



P: Temperature controller M1: Compressor M2: Fan

M2: Fan XP: Power Plug KA: Start relay C1: Starting capacitor

1. Installation and Start-up

1.1 Placement

1. For safe and operational reasons, the freezer must not be installed outdoors.
2. The freezer should be placed on a stable, level surface in a dry, well-ventilated room (max. 75% relative air humidity).
3. Never place the freezer close to heat sources such as cookers or radiators, and avoid placing it in direct sunlight.
4. Standard EN378 specifies that the room in which you install your freezer must have a volume of 1m³ per 8 g of hydrocarbon refrigerant used in the products, so as to avoid the formation of inflammable gas/air mixtures in the room where the freezer is located in the event of a leakage in the refrigeration circuit.

1.2 Room Temperature

The freezer is designed to operate at limited ambient temperatures according to the climatic category specified below.

- **Climate class optimum room temperature:**

N : +16 °C to +32 °C

1.3 Ventilation

It is important that the freezer is installed at a ventilated place and fresh air can circulate unhindered above, below and around it. The freezer requires at least 20cm free space at all around it.

2. Electrical Connection

Data regarding voltage, power and current are given on the nameplate which is at the back of the freezer.

Electrical connection must be made in accordance with the local regulations.

Warning!

- **This freezer must be earthed.**
- **If the socket and the plug are not the same type, please have the plug or socket replaced by a qualified electrician.**
- **Do under no circumstances use extension lead or adaptor.**
- **If the power cord is damaged, it must be replaced by the manufacturer, or service provider, or a qualified electrician in order to avoid a hazard.**

3. Temperature Controller

The freezer is equipped with the electronic controller.



3.1 Description

This KELD model is designed for many cooling applications. It has an input for temperature probe type PT. The probe temperature is displayed on the bright 3-digit display. The user is able to program 13 different parameters including shutdown temperature. The unit features error warning and password protection. Select between red, green or blue display color, and 115Vac, 230Vac, 24Vac/dc or 12Vac/dc power supplies.

3.2 List of Parameters

Code	Description	Units	Range	Default
SP	Shutdown temperature	Degrees	r1 to r2	-65
r0	Temperature difference	Degrees	0.1 to 20.0	3
r1	Minimum value for SP	Degrees	-200.0 to r2	-120
r2	Maximum value for SP	Degrees	r1 to 50.0	-25
P1	Temperature adjustment	Degrees	-20.0 to 20.0	0.0
A0	Alarm canceled	Degrees	0 to 99.9	2.0
A1	Alarm value	Degrees	0 to 99.9	0.0
A2	alarm disabled	Minutes	0 to 250	0
A8	Alarm delay	Minutes	0 to 250	0
H5	Password	Numeric	0 to 999	0
c0	Minimum downtime	Minutes	0 to 59	0
c2	Error status	Option	on or off	on
P4	Decimal display	Option	yes or no	Yes

3.3 Parameter programming

Set Point (SP) is the only parameter the user can access without code protection.

1. Press SET. SP text will appear on the display.
2. Press SET again. The real value is shown on the display.
3. The value can be modified with the UP and DOWN arrows.

4. Press SET to enter any new values.
5. Press SET and DOWN at the same time to quit programming or wait one minute and the display will automatically exit programming mode.

3.4 Access to all code protected parameters

1. Press SET for 8 seconds. The access code value 0 is shown on the display (unit comes with code set at 0 from factory).
2. With the UP and DOWN arrows, code can be set to user needs.
3. Press SET to enter the code. If code correct, the first parameter label is shown on the display (SP).
4. Move to the desired parameter with the UP and DOWN Keys.
5. Press SET to view the value on the display.
6. The value can be modified with the UP and DOWN arrows.
7. Press SET to enter the value and exit.
8. Repeat until all necessary parameters are modified.
9. Press SET and DOWN at the same time to quit programming or wait one minute and the display will automatically exit programming mode.

* The keyboard code can be reset to ZERO by turning off the controller and turning it on again while keeping the SET pressed.

3.5 Default Working

In case of probe error, the control performs start or stop according to C2.

3.6 Led Indication and Display Messages

The led OUT indicates if the load is connected or not. If the control is waiting the Minimum stopping time c0 to start a cool cycle the led flashes (50% ON 50% OFF).

- **Ooo**: Open Probe Error
- **---** : Short Circuit Probe Error
- **ALH**: The high temperature alarm
- **ALL**: The low temperature alarm

4. Defrosting

Items stored in the freezer and humidity of the surrounding air might gradually form frost and ice inside the freezer. A thin layer of the frost or ice will not affect the freezer performance,

but if the frost and ice are getting more and more, the energy consumption will get higher. Small amount of loose frost or ice can be removed by a plastic shovel.

Never use electrical or sharp tools as these may damage the inside chamber or shelves of the freezer.

When the frost or ice is accumulated too much, the freezer should be completely defrosted. It is recommended that defrosting is carried out after all stored items are moved out of the freezer.

Defrosting process as follows:

1. Power off the freezer, and leave all the freezer doors open to let the ice melt.
2. Use a soft dry cloth to wipe off the water inside of the freezer
3. Ensure the freezer inside clean and dry before powering on again.
4. Move items back to the freezer when the set temperature has been reached.

5. Cleaning and Maintenance

5.1 Cleaning

1. Before cleaning, cut off the power supply and remove all the stored items out of the freezer
2. The freezer is best cleaned by a clean, soft and dry cloth using mild water, if needed, neutral detergent can be used to clean the freezer. Never to scour the whole unit for cleaning. After cleaning, make sure the freezer get dried before powering on again. And please be careful to prevent water entering the control panel.
3. The sealing gaskets are recommended to clean regularly to prevent degradation and extend its span life. The sealing gaskets shall be cleaned by soft cloth with clean water.
4. If possible, it is also recommended that the tubular condenser in the refrigeration cabin to be cleaned in a regular basis to prevent too much dust and dirt affecting the refrigeration performance.
5. The ventilation grill also needs to be away from dust and dirt, and it is recommended to be cleaned with a vacuum cleaner in a regular basis.
6. If the freezer is not used for a long time, please switch off, empty and clean the unit before storage.

5.2 Energy Saving

1. Keep the inner doors closed as always
2. Make sure the inner doors and outer door closed properly after opening the freezer
3. Position the freezer in a cooler (min. 5°C), well-ventilated room.
4. Defrost the freezer at a regular time to prevent ice accumulating too much

6. Trouble Shooting

Fault	Possible Causes	Remedy
No control lamps are lit.	<ol style="list-style-type: none">1. The freezer is switched off.2. Power failure;3. the fuse is burned;4. the freezer is not plugged correctly.	<ol style="list-style-type: none">1. Check the plug and the socket.2. Check the power is connected correctly.3. Reset the fuse.
Temperature in freezer too high.	<ol style="list-style-type: none">1. The ventilation grill is blocked by dust and dirt.2. The inner doors and outer door are not closed properly.3. Too much items are placed into the freezer at a time.4. Too much frost and ice.5. The thermostat is set as a high point.	<ol style="list-style-type: none">1. Clean the grill to ensure unhindered air circulation.2. Close the inner doors and outer door properly and wait for a few hours.3. Wait a few hours to see if the temperature can reach the set point4. Defrost the freezer.5. Set the thermostat to a lower point.
Temperature in freezer too low	<ol style="list-style-type: none">1. The thermostat is set as a low point.2. The inner doors are prevented from closing tightly by items stored in the freezer.	<ol style="list-style-type: none">1. Set the thermostat to a higher point.2. Reposition items inside the freezer and make sure the inner doors are closed tightly.
The doors cannot be opened.	Vacuum pressure in freezer.	Wait 2-3 minutes until the vacuum pressure has been equalized and try again.
Vibration and unusual noise.	<ol style="list-style-type: none">1. The freezer is not placed on a stable and level surface.2. The freezer is placed too close to other equipment3. Fan is not run in a proper way	<ol style="list-style-type: none">1. Position the freezer on a flat, level and stable surface.2. Move other equipment away from the freezer3. Open the cabin and check if the fan is lack of lubricating oil

7. Safeguarding on the Environment

7.1 Packing

The packaging material is entirely recyclable, and marked with recycling symbol. Please follow the local regulations for scrapping.

7.2 Product

The freezer is made of recyclable material, while the refrigeration system is charged by flammable gases. Disposal of the entire freezer should be in accordance with the local environmental regulations and relative standards. Make sure the product is disposed by a qualified authority for recycling both freezer bodies and refrigerants.

7.3 Disposal

- Disposal of the freezer must comply with the local environmental regulations. Please contact your local authorities who shall guide you how to dispose or scrap it.
- This product is marked in accordance with EU Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).
- Ensuring the product is disposed in a correct way helps prevent potential negative impacts on the environment and human health.
- The symbol on the product or its accompanying documents means the product must not be disposed together with household refuse. Instead, it must be disposed at a proper recycling station.
- The product must be disposed in accordance with the current environmental regulations of the disposal of waste.
- For further information on the handling, recycling and reuse of this product, please contact your local authorities, qualified service providers, or the dealer from whom you bought the product.

Certificate of Quality

Checker: